

*Always one
step ahead!*



**What's new
in 2017**



Perio-Protection-System

Optimize
your practice
management

Active Oxygen Gel and Mouthwash for the treatment of Periodontitis and Peri-implantitis

OXYSAFE Gel

- Active oxygen gel
- Fights harmful bacteria and enhances wound healing
- Reduction of pocket depth
- Free of antibiotics and free of CHX
- Easy, precise and ready-to-use

OXYSAFE Liquid

- Active oxygen mouthwash
- Fights harmful bacteria and enhances wound healing
- Reduction of pocket depth
- Free of antibiotics and free of CHX

Article	REF
OXYSAFE Intro Kit, 3 x 1 ml gel, 3 x 250 ml Liquid	155 040
OXYSAFE Gel, 3 syringes of 1 ml	155 041
OXYSAFE Liquid, 250 ml	155 042

Miraject Flexi Tip

Application tips with Luer-Lock

- Blunt plastic needle pre-bend tip
- Very flexible and form stable
- For exact application and dosing of highly fluid pastes and liquids
- Single use only

Article	REF
Miraject Flexi Tip Intro Kit (3 Syringes 21G à 1 ml and 10 Single use tip)	254 240
Miraject Flexi Tip 21G, ø 0.6 – 0.8 mm x 24 mm (25 Single use tip), size thick	254 241
Miraject Flexi Tip 27G, ø 0.2 – 0.4 mm x 20 mm (25 Single use tip), size fine	254 242



Miraject Flexi Tip 21G ø 0.6 – 0.8 mm x 24 mm
Miraject Flexi Tip 27G ø 0.2 – 0.4 mm x 20 mm

Miratray 3 in 1

3 work steps in one - impression, counter impressionbite registration

- Time saving: With only one impression a double impression of upper and lower jaw is made
- Material savings: partial impression with direct correlation of the occlusion
- Patient friendly: faster and less gag reflex
- Sheer, tearproof net material
- Dimensionally stable with strengthened handle
- Color-coded

Article	Content	REF
Posterior, green	36 pcs.	101 270
Quadrant, yellow	36 pcs.	101 271
Anterior, purple	32 pcs.	101 272
Full Arch, blue	28 pcs.	101 273



Integra® Cart Smart

Cost-friendly 3 surface cart

- 2 universal electric sockets
- 2 trays with rim protection
- Maximum load per tray 10 kg
- Weight: 14.5 kg

measurements
35 x 78 x 45 cm

Article	REF
Integra Cart Smart	460 120



Back side



Our protective eyewear
Safe and stylish

iSpec Pure Fit II

Nonslip and comfortable

- Modern, sporty design
- Very comfortable and with soft temple and nose contact
- Optimal protection by twist side technology and nonslip
- Large field of vision due to duo-spheric glasses
- 27 g

Article	REF
blue glasses transparent	355 627
red glasses transparent	355 628



Protecto-Tip

Universal protection cap for instruments and inserts

- Caps made of flexible plastic
- Protects pouches from punctures
- Avoids risk of injuries according to RKI
- Simple and fast application
- For almost all forms of instruments
- Also recommended for Cavitron Inserts
- Autoclavable

Article	Content	REF
2.5 mm x 28 mm	10 pcs.	450 110
4.0 mm x 28 mm	10 pcs.	450 111



H&W OTG II

Over-the-glass panorama overspecs

- No contact of temples with the prescription glasses
- Good coverage of the eye area
- 38 g

Article	REF
H&W OTG II	355 618



H&W Panorama II

For slender faces

- Temple with ergonomic inclination
- 25 g

Article	REF
H&W Panorama II	355 617



HW Super Fit UV (OTG)

Over-the-glass UV Blocker

- Suitable for works with polymerization lamps (390-525 nm)
- Softpad technology assures ideal, ergonomic fit
- Unobstructed view to the sides
- Temple length adjustable
- Coating: Outside extremely scratch-resistant, inside permanently mist-free
- 38g

Article	REF
H&W OTG UV	355 626





*miradent -
for a bright and shiny smile*

Mirawhite® shine

Effective tooth whitening even for sensitive teeth

- Tooth whitening by up to five shades
- Pen for easy and precise application
- Tooth whitening with active-oxygen (Ardox-X® Technology), free of peroxide
- Gentle on enamel and gingiva
- Period application 7 – 14 days, twice a day

Article	REF
Mirawhite shine, 1.8 ml	630 194
Mirawhite shine, with 6 pcs. Display	635 212



XyliPOP®

Dental care lollipop sweetened with Xylitol – sugar free!



- With tooth-friendly Xylitol and calcium lactate
- Natural ingredients
- Free of aspartame, sorbitol, lactose and gluten
- Motivation and fun in childrens' prophylaxis
- A "dental care" sweet for children and adults alike
- Fruity-fresh strawberry flavor

Article	Content	REF
XyliPOP Lolli Erdbeere	6 g	630 208
XyliPOP Lolli Erdbeere	40 pcs. x 6 g	635 211



Funny Animals

Hygienic toothbrush holder

- Easy snapping in place and removing of toothbrush
- Patented snap lock
- Two suction cups for attaching to all even surfaces
- Ventilation slots ensure an ideal air exchange
- Fitting for any off-the-shelf toothbrush - for children aged 3 years and older

Article	REF
Funny Cow	630 099
Funny Giraffe NEW	630 220
Funny Panda NEW	630 222
Funny Lion	630 102
Funny Elephant	630 100
Funny Hippo NEW	630 221



I-Prox® L


Interdental brushes with L-shaped handle

- L-shaped handle and 100° angle facilitate handling
- Hygienic protective cap
- Ideal on the go
- Six sizes from xx-fine - large (Color-coded)

Design	Insert Cylinder	Cleaning Cylinder	REF
xx-fine	Ø 0.4 mm	Ø 1.8 mm	630 210
x-fine	Ø 0.5 mm	Ø 2.0 mm	630 211
fine	Ø 0.6 mm	Ø 2.5 mm	630 212
medium	Ø 0.7 mm	Ø 3.5 mm	630 213
conical	Ø 0.8 mm	Ø 2.5 - 5.0 mm	630 214
large	Ø 0.8 mm	Ø 5.0 mm	630 215



Best Price



Your partner
for dental specialities,
ideas and service.



Häger & Werken GmbH & Co. KG
Ackerstraße 1, 47269 Duisburg, Germany
Tel. +49 (203) 99269-0
Fax +49 (203) 299283
E-Mail info@hagerwerken.de
Internet www.hagerwerken.de

Please see our general terms and conditions at www.hagerwerken.de

