




Science.
Applied to Life.™

Raising the bar
with stability
and control.

3M™ Incognito™ Brackets with Tip Bar
Invisible, intelligent, individualized.





Incognito™ Bracket with Tip Bar

“ I evaluated the new Incognito™ Brackets with Tip Bar, with mesial and distal extensions to the slot on the lower canines and upper 3-to-3. I found that the slot extension helps to **control and correct the angulation** of those teeth in the early stage of development. This product improvement will be essential to my future cases, especially considering the **reduction in chair and treatment time.** ”

– Dr. Esfandiar Modjahedpour

Latest improvement in a series of ongoing innovations.

Creation of the 3M™ Incognito™ Brackets with Tip Bar continues to enhance lingual treatment with the 3M™ Incognito™ Appliance System, delivering targeted results through a new feature that improves angulation control. Extensions on the mesial and distal sides of the lingual bracket elongate the bracket slot to increase the stability of the interplay between the brackets and archwires.

Setting a new standard for completely aesthetic, comfortable and individually customized orthodontic options, the Incognito Bracket with Tip Bar design ...

- Enables utilization of the Incognito Bracket with Tip Bar for as many teeth as possible, even for crowded cases.¹
- Offers same or better patient comfort,¹ ease of use,¹ and handling¹ as the conventional Incognito Bracket.
- Is available on lower canines and upper 3-to-3.
- Provides even more angulation movement when combined with a 3M™ AlastiK™ Lingual Ligature.¹

REFERENCES

1. Data on File

Raise the bar today!

Visit 3M.com/Incognito
or contact your 3M
Representative.

Incognito™
Appliance System

3M, AlastiK, and Incognito are trademarks of 3M. Used under license in Canada.
© 3M 2017. All rights reserved. 70-2013-0838-7 1701

3M