




Lingual Orthodontic Products and others


- Orthodontic Specialities
- Special Lingual Pliers
- Pliers and Instruments
- Bonding Supplies
- Microetcher
- Wires, Ligatures, Elastics
- Photography, Accessoires
- Miscellaneous



SPECIAL LINGUAL PLIERS

Distal Cinch Back Plier – Long Handle Serrated jaws, works perfectly with round NiTi Wires

Item-No. 557-0005M



New Style, bends NiTi wires distal to the molar tube. Long handle (total length: 142 mm).

Working range, NiTi .012 – .017 x .025

Same purpose as 557-0005M, but without serration.

Working range, NiTi .012 – .017 x .025



Item-No. 557-0005

Three Jaw Plier (Aderer)

Item-No. TASK-60-313


For permanent bending
of all NiTi wires.

**Maximum bending
Capacity: Up to .028**



Lingual Utility Plier

Item-No. 557-0001



Slim, gracile plier, especially suited for the molar; lingual, but also vestibular. Total length: 145 mm.

SPECIAL LINGUAL PLIERS

Lingual Utility Plier Long Handle

Item-No. 557-0001XL

Long handle (total length: 175 mm).




Lingual Debonding Plier – Modified

Item-No. 553-0431P

Lingual Debonding Plier.
Ideal for debonding of lingual brackets in the front area.

Note: Now with lateral inbus screw to fix or replace the pin. It is very easy to fix or replace the pin.



Replacement Pin

Item-No. 553-0430

Debonding Plier for Molars

Item-No. 557-0006A



Debonding Plier. Modified special cutting edge to debond lingual brackets and tubes in the posterior area. Attention, no cutter! This modified Debonding Plier isn't suitable for cutting wires. The black colored handle prevents from confusion with cutters.



SPECIAL PLIERS

Ezcurra Lingual Mini Distal End Cutter with Hold – Slim Long Handle

Item-No. 550-0510

NEW!

Replaces the plier
AEZ-803-0510.




One Piece Mini Distal End Cutter with Hold. These cutter has the additional feature of a mini tip designed for easier access distal to the buccal tube. The blade design enables the clinician to make cuts within 0.5 mm of the buccal tube, while the „hold“ feature prevents the cut-off portion of the wire from springing loose. Smallest head of Distal End Cutters.

Maximum cutting capacity: Wires up to .019 × .025

Distal End Cutter Safety Hold – Flush Cut

Item-No. TASK-30-552





Lingual Distal End Cutter

Item-No. HF678-700

Occlusal view on the Lingual Distal End Cutter cutting the excess archwire distal to the buccal tube.



Cutting capacity: .012 to .018



SPECIAL LINGUAL PLIERS

Reversed Ligature Cutter – 40° Angle

Item-No. 550-0410




Ligature Cutter, same model as 550-0400 (see on catalog page B2 above), but with internal cutting surface. 40° angle. Cuts the ligature wire distal to the bracket.

Maximum cutting capacity: Ligature wires up to .008

Pin and Ligature Cutter

Item-No. 550-0605S




Slim Pin and Ligature Cutter for easy lingual access. Cuts the ligature wire mesial to the bracket.

Maximum cutting capacity: Ligature wires up to .015

Micro Mini Pin and Ligature Cutter – Long Handle

Item-No. 550-0107L




Long handle (total length: 135 mm).

Maximum cutting capacity: Ligature wires up to .012

SPECIAL PLIERS

Intraoral Detailing Plier – Left

Item-No. 551-0643A




Ideal for finishing bends, bypasses and 2nd order bends.

Makes 0.75 mm bends (in/out, up/down).

Wires up to .018 x .025

Intraoral Detailing Plier – Right

Item-No. 551-0644A



Ideal for finishing bends, bypasses and 2nd order bends.

Makes 0.75 mm bends (in/out, up/down).


Wires up to .018 x .025



Weingart Utility Plier – Hardened Steel Insert

Serrated inserts made from hardened steel.


Item-No. TASK-60-501T



Weingart Utility Plier – Hardened Steel Insert

Precise closing, fine tips, very sturdy. With hardened steel inserts.


Item-No. 550-6556



Weingart Utility Plier Fine – Hardened Steel Insert

Fine tips. Serrated inserts made from hardened steel.

Item-No. 550-0502



Weingart Utility Plier XL – Long Handle

Weingart Utility Plier with fine, serrated tips from special hardened steel. Long handle (total length: 155 mm) design allows easy posterior access.


Item-No. 550-0505XL



Weingart Plier Gracile – Tungsten Carbide Insert

Gracile tips. Serrated inserts made from tungsten carbide.

Item-No. 550-0503A



Pin and Ligature Cutter Miniature – Hardened Steel Insert

Maximum cutting capacity: Soft wires up to .015


Item-No. TASK-60-605



TASK Ligature Cutter – Hardened Steel Insert

Ligature Cutter, same model as 550-0410 (see page A5), but with external cutting surface. 40° angle.
Maximum cutting capacity:
Ligature wires up to .008


Item-No. 550-0400



Hard Wire Cutter – Hardened Steel Insert

Maximum cutting capacity:
Hard wires up to .020, round

Item-No. TASK-60-604H



Distal End Safety Cutter Slim – Tungsten Carbide Insert

Long handle.
Maximum cutting capacity:
Rectangular wires up to .022 x .028
Round wires up to .020

Item-No. TASK-30-550SL



Lingual Distal End Cutter – Flush Cut

With steel inserts. Long handle, easier access to lingual and posterior area. Cuts without rotational force. Safely cuts & holds distal end of wire. Cutter can be sterilized after removing rubber ring. Ring to be sterilized separately.

Maximum cutting capacity:
Rectangular up to .019 x .025
Round up to .020

Item-No. 550-0550FL



Replacement Rubber Retaining Rings: 10 pieces per pack

Service life of the ring can be extended by rotating it after use. Ring to be sterilized separately.

Item-No. 550-0550-FL/R



Distal End Cutter Long Handle – Tungsten Carbide Insert

With tungsten carbide inserts. Length: 150 mm.
Wires from .014 up to .021 x .025.


Item-No. 550-632/9TC



Distal End Safety Cutter Slight Head – Tungsten Carbide Insert

With tungsten carbide inserts.
Wires from .012 – .021 x .025.

Item-No. 550-528L



Lingual Distal End Cutter – Flush Cut

Long handle, easier access to
lingual and posterior area.
Cuts without rotational force.
Safely cuts & holds distal end of wire.
Maximum cutting capacity:
Up to .019 x .025

Item-No. 557-0004A



Direct Bond Attachment Removing Plier – Hardened Steel Insert

With hardened steel inserts.

Straight Item-No. TASK-TS-25
Angled Item-No. TASK-TS-26



Direct Bond Attachment Removing Plier Angled – Hardened Steel Insert

With hardened steel inserts.


Item-No. TASK-60-334T



Tweed Arch Bending Plier – Hardened Steel Insert

With Inserts. Maximum bending capacity:
Rectangular wires up to .022 x .028

Item-No. TASK-60-303T



Tweed Arch Bending Plier

Wires up to .022 x .028


Item-No. TASK-60-303



Loop Forming Plier

Maximum bending capacity:
Wires up to .022

Item-No. TASK-60-305



Omega Loop Forming Plier

Wires up to .020


Item-No. TASK-60-304S



Kim Plier with Cutter

Wires up to .016 x .022

Item-No. TASK-60-430CK



Hook Crimping Plier – Angled




Item-No. TASK-60-446



Hook Crimping Plier – Straight



Item-No. TASK-P220



Heavy Three Jaw Plier (Aderer)

Wires up to .036

Item-No. TASK-60-313L



Universal Distal End Cutter


With normal or long handle. With hold.

Long handle, wire size: .021 x .025

Normal handle, wire size: .021 x .025

Item-No. HF678-101L

Item-No. HF678-101



Micro Mini Ligature Cutter – 15° Angle

Wire size up to .012


Item-No. HF678-109




Lingual Pin and Ligature Cutter

Wire size up to .012

Item-No. HF678-701



Occlusal view to the Lingual Pin and Ligature Cutter cutting the excess steel ligature.




Lingual Weingart Plier

90° bend, up to .022

Item-No. HF678-702-90



Occlusal view to the Lingual Weingart Plier guiding the archwire into the bracket slot.




Lingual Mathieu Plier

Item-No. HF678-705



Occlusal view to the Lingual Mathieu Plier tying steel ligature to hold the arch wire into bracket slot.



Lingual Bracket Removing Plier

Item-No. HF678-703



Band Pusher Scaler Small

Item-No. 300-23-003



Micro Pin and Ligature Cutter

Narrow and slim head designed for easy access to hard-to-reach areas and for cutting soft wires, pins and elastics. Corrosion resistant metal inserts.

Soft wires up to .015"

Item-No. 300-5002



Pin and Ligature Cutter

Straight head, standard handle, abrasion and corrosion resistant metal inserts.

Soft wires up to .015"

Item-No. 300-5006



Ligature Cutter 40°

Standard handle, abrasion and corrosion resistant metal inserts, Metal inserted jaws ensure abrasion constant performance and reliability of cutting edges.

Soft wires up to .015"

Item-No. 300-5009



Hard Wire Cutter

Straight head, standard handle, abrasion and corrosion resistant metal inserts.

Wires up to .022 x .028"

Item-No. 300-5018



Distal End Cutter With Safety Hold – Long Handle

Long handle, tungsten carbide inserted jaws.

Minimum: Wires up to .014

Maximum: Wires up to .021 x .025"

Item-No. 300-5028



Distal End Cutter With Short Head

Narrow and slim head,
tungsten carbide inserted jaws.

Minimum: Wires up to .012

Maximum: Wires up to .022 x .028"

Item-No. 300-5029



Flush Cut Distal End Cutter Safety Hold

Narrow and slim head, standard handle,
tungsten carbide inserted jaws.

Minimum: Wires up to .014

Maximum: Wires up to .021 x .025"

Item-No. 300-5026




Flush Cut Distal End Cutter Safety Hold – Long Handle

Long handle, tungsten carbide inserted jaws.

Minimum: Wires up to .014

Maximum: Wires up to .021 x .025"

Item-No. 300-5030



Weingart Plier Slim

Standard handle. Slim jaws provide easy access
to difficult areas. Bends and holds round,
square and rectangular wires up to .030"

Item-No. 300-5032



Weingart Plier Angled Short Beak

Short beak, standard or long handle, precise
and accurate work. Narrow and slim head
allowing easy access to all areas.

Wires up to .030"

Standard Handle

Long Handle

Item-No. 300-5040

Item-No. 300-5040L

Complete Lions pliers catalog on request.

SCHWERT PLIERS



SCHWERT
SINCE 1874
1918 - 2018



Ligature Cutter

Length: 11.5 cm.
With HCS inserts and screw lock.
Soft wires up to .016/0.4 mm

Item-No. 550-5248



Ligature Cutter

Length: 12.0 cm.
With HCS inserts and screw lock.
Soft wires up to .020/0.5 mm


Item-No. 550-5253



Wire Cutter

Length: 12.0 cm.
With HCS inserts and screw lock.
Hard wires up to .028/0.7 mm

Item-No. 550-5254



Mini Distal End Cutter

Length: 12.0 cm.
With extra hard steel inserts.
Hard wires up to .022 x .028/0.56 mm x 0.70 mm

Item-No. 550-5230


SCHWERT LABOR PLIERS



Flat Nose Plier Crovana

Length: 12.5 cm. Serrated tips.
Wires up to .028/0.7 mm


Item-No. 550-5022



Lingual Tube and Bracket Tweezer Fine

Angulated fine, self-locking. Slim tips for positive grip, one tip fits into the slot of any lingual or buccal tube.


Item-No. 553-0503



Lingual Tube and Bracket Tweezer Offset Fine – Modified

The same purpose as above, but with offset.

Item-No. 553-0504



Ligature Director (Tucker) Double Ended

Slot in tips is chamfered to avoid damaging elastics. Both ends are useful for intraoral finishing bends.

Item-No. 560-0010



Diamond Ligature Director

The diamond coated tip prevents slipping and simplifies the correct placing of ligatures.

Item-No. TASK-74822



Probe Short Explorer


5 mm tip. Slightly serrated handle.

A: thin handle, \varnothing 4.5 mm

B: thick handle, \varnothing 6.0 mm

Item-No. 553-3007


Item-No. 553-3008



Angled Explorer Single Sided

Thicker handle \varnothing 8,1 mm.

Item-No. 300-33-012



Spatula Heidemann

Spatula, elastic, extra flexible. 160 mm

Item-No. ADE-420



Mosquito Forceps Straight


For easy placement of elastic modules, serrated tips for positive gripping.

Item-No. 554-1020

Mosquito Forceps Curved

For easy placement of elastic modules, serrated tips for positive gripping.

Item-No. 554-1021



Smaha Forceps Needle Holder Straight – Mathieu Style

Allows access to lingual and posterior areas. Very slim beak, with serration.


Item-No. 554-1026



Smaha Forceps Needle Holder Curved – Mathieu Style

Allows access to lingual and posterior areas. Very slim beak, with serration.


Item-No. 554-1025



Elastic Needle Holder Hook Tip – Without Serrated Tips

With small hook on the tip. For easy placement of elastic modules, without serrated tips.


Item-No. 550-0216



Needle Holder Lobster – Mathieu Style

Mathieu style ligating plier.
Excellent for placement of elastic ligatures.
Locking handle provides non slip grip.



Item-No. OL-200-701MK



Needle Holder – Mathieu Style

Standard, narrow serrated tips.



Item-No. OL-200-504PS



Needle Holder – Mathieu Style

Standard, medium serrated tips.



Item-No. 580-0026A



Needle Holder – Mathieu Style

Standard, medium serrated tips,
with TC inserts.

Item-No. OL-200-503PB



Needle Holder – Mathieu Style

Standard, wide serrated tips.

Item-No. 580-0026

PLIER RACKS



Pliers Carousel

Designed for easy reach and visibility.
The inner clear tube holds additional straight instruments. Clear plastic.

127 mm diameter footing, 38 mm inner ring
102 mm diameter outer ring, 140 mm high

1 piece

Item-No. OL-200-240



Pliers Rack

1 piece

Item-No. HF-200



Upright Pliers Rack

Upright Pliers Rack with cover holds 15 pliers.
White, horizontal storage tray at front.

L 190 mm x W 152 mm x H 164

1 piece

Item-No. OL-200-243

INSTRUMENT CARE



ILS Lubricant Spray

Hu-Friedy Care spray for free-standing hinge instruments. To avoid corrosion. Application after each cleaning.

236 ml bottle

Item-HF-ILS

IPS Penetrating Oil

Hu-Friedy penetration oil for hinged instruments. Application after each cleaning. Rust is not removed. Agility can be restored. Corroded instruments must be removed from the processing circuit.

236 ml bottle

Item-No. HF-IPS



IMS Shine reNEW Instrument Wipes

Effectively removes stain and rust, restoring the luster to instruments. Keeps hinged instruments moving freely by dissolving built-up residue. Eliminates the need for abrasive scrubs. Quick and simple to use.

20 wipes

Item-No. HF-IMS-1455



Instrument Oil-Care Pen

Precision point oilers for lubrication and conservation of pliers hinges. Clean and precise. The synthetic food grade, high-tech oil is transparent and odourless. In contrast to oils based on mineral oils, it is toxicologically harmless. All the ingredients are USDA approved and comply with the FDA or the German pharmacopoeia (DAB). It also provides outstanding high-temperature stability (flash point 260° C), excellent wear protection as well as lubrication and penetrating properties. The oil is suitable for all sterilisation procedures. During sterilization, the oil spreads evenly in the joint and rinses out possible dirt residues. No oiling is necessary. Drops are applied optimally at room temperature.

2 pieces in case

Item-No. OL-200-990



Plastic Conditioner

For bonding to composite restorations and acrylic appliances.

10 cc bottle

Item-No. 740-3501



Porcelain Conditioner

For bonding to ceramic crowns.

8 cc bottle

Item-No. 740-3502



Metal Primer

For bonding to gold, stainless steel and amalgam.

8 cc bottle

Item-No. 740-3503



Porc Etch with Barrier Gel

Porcelain Etching Agent with Barrier Gel.

**Porc Etch: 7 g bottle
with Barrier Gel: 3 cc bottle**

Item-No. 740-3504

Barrier Gel: 3 cc bottle

Item-No. 740-3504/B



Maximum Cure Sealant A

For indirect bonding with silicon tray. Chemically curing. Store in the fridge.

20 g bottle

Item-No. 740-3507

7 g bottle

Item-No. 740-3505



Maximum Cure Sealant B

For indirect bonding with silicon tray. Chemically curing. Store in the fridge.

20 g bottle

Item-No. 740-3508

7 g bottle

Item-No. 740-3506

We carry all Reliance Products!



Maximum Cure Cooler

Very efficient for cooling and storage of Maximum Cure A & B at the chair side. For two 20 g bottles. Blue, special aluminium alloy.

Item-No. 740-3600



Nexus III

Chemical and light cure composite. Clear. Double syringe 5 g, including 10 mixing tips. Recommended for rebonding of lingual brackets.

5 g syringe, 10 mixing tips

Item-No. 740-4336N

Nexus III Automix Tips

Refill Automix Tips.

50 pieces per pack

Item-No. 740-4338N



Excite F DSC

Refill Regular. Soft Touch™ single dose. For bonding with clear transfer trays. Contains fluoride.

20 x 0.1 g per pack

Item-No. 740-1262

50 x 0.1 g per pack

Item-No. 740-1261



Rely-a-Bond

Single-phase orthodontic bonding system. For bonding to metal, composite and ceramic brackets.

Paste: 4 x 3.5 g syringes

Item-No. 740-3001

Paste: 14 g jar

Item-No. 740-3002

Primer: 16 cc bottle

Item-No. 740-1002



Phase II – Paste A / Paste B

3 x 6 g in syringes Paste A

Item-No. 740-4001

3 x 6 g in syringes Paste B

Item-No. 740-4002



Phase II – Paste A

18 g in jar

Item-No. 740-4005

Phase II – Resin A

7 g bottle

Item-No. 740-4003

20 g bottle

Item-No. 740-4007



Phase II – Paste B

18 g in jar

Item-No. 740-4006

Phase II – Resin B

7 g bottle

Item-No. 740-4004

20 g bottle

Item-No. 740-4008



Light Bond and Light Bond Thin

Curing recommended with L.E.D. lamps.

Wavelength: 440 – 480 nm.

Kit: 4 x 5 g Light Bond Paste in push syringe;

2 bottles Sealant Resin à 3 cc;

23 cc Liquid Etchant

Item-No. 740-7004

**1 x 5 g syringe Light Bond Paste Thin
for rebonding**

Item-No. 740-7006/TH

4 x 5 g syringes Light Bond Paste

Item-No. 740-7001

3 cc bottle Light Bond Sealant

Item-No. 740-7002



Ultra Band-Lok

Light activated single paste. Band cement.

Ivory: 1 x 5 g syringe

Item-No. 740-5008

Blue: 1 x 5 g syringe

Item-No. 740-5007



Pro Seal for Halogen Lamp

Light cure fluorescing agent. Highly filled sealant. Contains fluoride. Halogen or high intensity curing lamps: wavelength 390 – 440 nm.

6 cc bottle

Item-No. 740-7020



Pro Seal L.E.D.

Light cure filled sealant. Contains fluoride.

L.E.D. curing lamps: wavelength 440 – 480 nm.

6 cc bottle

Item-No. 740-7025



Assure / Assure Plus

Universal bond enhancer for wet or dry environment.

Assure 6 cc bottle

Item-No. 740-6501

Assure Plus 6 cc bottle

Item-No. 740-6503



Light Cure Retainer Paste

Flowable, light cure orthodontic adhesive.

5 cc syringe

Item-No. 740-7011



Flow Tain LC / Flow Tain LC, LV

Flowable, light-cured composite. Ideal for lingual retainers, good application by needle tip.

Flow Tain LC:

1,5 g syringe + 10 tips

Item-No. 740-7000

Flow Tain LC, low viscosity:

1,5 g syringe + 10 tips

Item-No. 740-7000/LV



Etch Rite

Viscous phosphorus gel (38%) for the etching of enamel and dentin. Effects easy application to multi-banded appliances.

4 syringes, à 1.2 ml

Item-No. OL-002-2100

24 syringes, à 1.2 ml

Item-No. OL-002-2130

Jumbo-Kit, 2 x 25 ml (64 g)

5 x 3 ml empty syringes, 50 tips

Item-No. OL-002-2140

Jumbo-Refill, 2 x 25 ml (64 g)

Item-No. OL-002-2150

50 Needle tips, curved

Item-No. OL-002-2170



Etchant Agent Liquid


37% phosphoric acid. For use on enamel to prepare the surface for bonding.

9 g bottle

Item-No. 740-4500

18 g bottle

Item-No. 740-4501



Dimensions:
L: 265 mm
Ø: 27 mm
Weight: 205 g

Curing Light

Model: C01-D (Dual Range)
Fiber Optic Light Guides: 11 mm diameter
Wavelength: 390-430 nm / 440-480 nm

Curing Light Black	Item-No. C01-DBK
Curing Light Light-Blue	Item-No. C01-DLB
Curing Light Light-Green	Item-No. C01-DLG
Curing Light Light-Lavender	Item-No. C01-DLLV
Curing Light Lime	Item-No. C01-DLM
Curing Light Lavender	Item-No. C01-DLV
Curing Light Orange	Item-No. C01-DOR
Curing Light Pink	Item-No. C01-DPK
Curing Light White	Item-No. C01-DWT
Curing Light Yellow	Item-No. C01-DYL



Rely X Unicem 2 – Automix

- Clinically tested self-adhesive cement, with a very low risk of postoperative sensitivities.
- Unique self-adhesive chemistry for very high bond strength, eliminating the need for etching, priming and bonding.
- Limited discolouration, ensuring long-lasting aesthetics.

**8,5 g base paste/catalyst paste,
15 automix tips + 5 endo tips**

Item-No. PP-56846



Concise

Bonding Kit.

**18 g paste A, 18 g paste B,
6 cc resin A, 6 cc resin B,
9 cc etching liquid,
200 applicator sponges,
2 blocks, 50 mixing spatulas,
3 application brushes**

Item-No. PP-196-001




Transbond IDB

Transbond IDB Pre-Mix Chemical Cure Adhesive.

**10 g component A + B,
200 applicators,
20 disposable containers**

Item-No. PP-712-120



Mini Mold Button Bonder Kit

Creates a functional lingual button without impression by using light cure composite.

10 tips + handle per pack

Item-No. 625-0701



Mini Mold Habit Control Spike Bonder Kit

Forms a small conical spike on the lingual of the incisors. Great for correcting thumb sucking and tongue thrusting with minimal chair time.

10 tips + handle per pack

Item-No. 625-0704



Mini Mold Bite Ramp Kit

Creates a bite-opening ramp on the lingual of the upper anterior teeth by using light cure composite.

8 tips (3 mm) + handle per pack

Item-No. 625-0702

8 tips (5 mm) + handle per pack

Item-No. 625-0702/5



Mini Mold Starter Kit

Creates 9 different shapes of composite attachments, using light cure composite. Wire Bonder (large and small), Bracket Bonder, Button Bonder, Habit Control, Molar Pad, Bite Ramp (3 mm and 5 mm), Tube Bonder.

9 different molds (3 ea.) + 2 handles

Item-No. 625-0700



Mini Mold Handle

1 handle per pack


Item-No. 625-0709



Monoart Mixing Cups White

5 x 50 pieces per pack

Item-No. 735-1041



Monoart Flexible Applicators Blue

2 x 100 pieces per pack

Item-No. 735-1025



Parotisroll

Size 1: 80 mm, \varnothing 9 mm

100 pieces per pack

Item-No. OL-210-0001

Size 3: 100 mm, \varnothing 10 mm

100 pieces per pack

Item-No. OL-210-0003



Reliance Brushes Short

250 pieces per pack


Item-No. 740-4025



Reliance Brush Handle

1 piece per pack

Item-No. 740-4026



Reliance Brushes Long

250 pieces per pack

Item-No. 740-4250

DANVILLE MICROETCHER THE ORIGINAL



Before Microetcher II
Foil Mesh Bracket Base



After Microetcher II
Foil Mesh Bracket Base



Microetcher IIA
Completely autoclavable



Microetcher II

**Patient Face-Mask
see on page F4!**



Microetcher™ Model II

Only the nozzle is autoclavable. 360° rotatable and interchangeable nozzle. Only air pressure needed.

Approx. 6 bar using compressed air. Item-No. 850-9001



Microetcher™ Model IIA

High-grade steel, Intraoral Sandblaster. Fully autoclavable and 360° rotatable, interchangeable nozzle, permanently affixed jar lid. Only air pressure needed.

Approx. 6 bar using compressed air. Item-No. 850-9003



Microetcher™ Model CD

Allows to plug the unit directly into the KaVo or Sirona handpiece connector.

Pressure output: 3 bar.


**Sirona
KaVo
Bienair**

**Item-No. 850-9001CD-S
Item-No. 850-9001CD-K
Item-No. 850-9001CD-B**



Replacement Jar for Microetcher

Microetcher II: Jar with rubber seal in the lid	Item-No. 850-9004
Microetcher IIA: Replacement jar	Item-No. 850-9005




Quick Disconnect Kit

Danville Quick Disconnect Kit for Microetcher. For connection to compressed air supply. Full pressure air line is cut and slid onto T-Connector.

Quick Disconnect Kit	Item-No. 850-9002
Quick Disconnect Kit without male part	Item-No. 850-9007

Using the Quick Disconnect Kit





Aluminium Oxide

Air abrasion powder for sandblasting.

White: 50 micron (intraoral)

1000 g bottle

Item-No. 850-9006

Grey: 90 micron for laboratory

454 g bottle (extraoral)

Item-No. 850-9006/1



CoJet™ Sand

Blast-coating agent 30 micron, applied with sandblaster (Microetcher), effects a stronger bonding with metal, ceramic and composite surfaces.

40 g bottle

Item-No. 850-0962

Espe™ Sil

Directly after coating, dose Espe Sil in a dappen dish and use a clean brush to wet the entire surface with the solution.

8 cc bottle

Item-No. 850-0803




MacrocabPlus Dust Cabinet, Danville

Large self-contained dust cabinet for the lab, 220 V.

Size in mm: 201 h x 318 w x 267 d

Item-No. 850-9000P



MicrocabPlus Dust Cabinet, Danville

Small self-contained dust cabinet for the orthodontist's chair, 220 V.

Size in mm: 201 h x 221 w x 267 d


Item-No. 850-9010



Shield for MicrocabPlus

10 pieces

Item-No. 850-9034



Filter for MicrocabPlus

Item-No. 850-9032



Wall Mounting Kit for MicrocabPlus

Optional for wall mounting.
Microcab easily removable for cleaning.

Item-No. 850-9010W



Quick-Disconnect Male Part for Microetcher

Connection for the Microetcher tube.

1 piece per pack

Item-No. 850-9009



Nozzle for Microetcher II and IIA

.048, 60°.

Complete replacement nozzle

Item-No. 850-9011

Nozzle for Microetcher IIA – Slim Profile

.048, 60°.

Complete replacement nozzle

Item-No. 850-9012



Replacement Carbide Tip for Microetcher II and IIA

Tungsten carbide replacement tip .048
for Microetcher II & IIA.

1 piece per pack

Item-No. 850-9008

4-Hole



Connectors for Microetcher II and IIA

4-Hole-Connector

Item-No. 850-9041

KaVo-Connector

Item-No. 850-9023

Sirona-Connector

Item-No. 850-9024



KaVo

Sirona



Air Compression Coupling for Microetcher II and IIA

Quick-Connect. 7.2 mm for laboratory.

1 piece per pack

Item-No. 850-9027



Rounded „Neck“



Bonding Adhesive Remover

Carbide bur, conical round, right angle. Special shape, with rounded blades to the „neck“. Removes residual adhesive without enamel damages.

Length of head: 4.4 mm

Ø .016 5 pieces per pack

Item-No. 820-9016



Bonding Adhesive Remover

Carbide bur, conical round, friction grip, turbine.

Length of head: 4.4 mm

Ø .016 5 pieces per pack

Item-No. 820-9016FG



Rounded „Neck“

Bonding Adhesive Remover

Carbide bur, conical round, friction grip, turbine.

With rounded blades to the „neck“.

Length of head: 4.7 mm

Ø .016 5 pieces per pack

Item-No. 820-9017FG



Lingual Bonding Adhesive Remover

Carbide bur, egg-shape, right angle, ideal bond remover for concave surfaces.

Length of head: 4.2 mm

Ø .023 5 pieces per pack

Item-No. 820-9379



Lingual Bonding Adhesive Remover

Friction grip, turbine.

Length of head: 4.2 mm

Ø .023 5 pieces per pack

Item-No. 820-9379FG



Bonding Adhesive Remover

Carbide bur, conical round.

Length of head: 6.0 mm

Ø .016 5 pieces per pack

Item-No. 820-9282



Bonding Adhesive Remover

Friction grip, turbine, conical round.

Length of head: 6.0 mm

Ø .016 5 pieces per pack

Item-No. 820-9282FG



Diamant Round

Right angle. Diamant bur, for grinding any type of metal attachments / materials.

Ø .035 **5 pieces per pack**

Item-No. 820-9835

Ø .023 **5 pieces per pack**

Item-No. 820-9823



Diamant Round

Friction grip, turbine.

Ø .035 **5 pieces per pack**

Item-No. 820-9835FG

Ø .023 **5 pieces per pack**

Item-No. 820-9823FG

LINGUAL WIRES

G4™ NiTi Lingual Arches

Archforms are designed and adapted precisely to lingual anatomy. Various alloys are offered to maximize treatment efficiencies in the narrow interbracket dimensions of lingual orthodontic treatment.

2 pieces per pack

Wire Size	Small Upper	Small Lower	Medium Upper	Medium Lower	Large Upper	Large Lower
.010	GH-G4LSU010	GH-G4LSL010	GH-G4LMU010	GH-G4LML010	GH-G4LLU010	GH-G4LLL010
.013	GH-G4LSU013	GH-G4LSL013	GH-G4LMU013	GH-G4LML013	GH-G4LLU013	GH-G4LLL013
.016	GH-G4LSU016	GH-G4LSL016	GH-G4LMU016	GH-G4LML016	GH-G4LLU016	GH-G4LLL016


M5™ NiTi Lingual Arches


Archforms are designed and adapted precisely to lingual anatomy. Various alloys are offered to maximize treatment efficiencies in the narrow interbracket dimensions of lingual orthodontic treatment.

2 pieces per pack

Wire Size	Small Upper	Small Lower	Medium Upper	Medium Lower	Large Upper	Large Lower
.016	GH-M5LSU016	GH-M5LSL016	GH-M5LMU016	GH-M5LML016	GH-M5LLU016	GH-M5LLL016
.016 x .022	GH-M5LSU1622	GH-M5LSL1622	GH-M5LMU1622	GH-M5LML1622	GH-M5LLU1622	GH-M5LLL1622
.017 x .017	GH-M5LSU1717	GH-M5LSL1717	GH-M5LMU1717	GH-M5LML1717	GH-M5LLU1717	GH-M5LLL1717
.017 x .025	GH-M5LSU1725	GH-M5LSL1725	GH-M5LMU1725	GH-M5LML1725	GH-M5LLU1725	GH-M5LLL1725

G4™ + M5™ Archforms






Crimpable Arch Stops (open)

.018 x .025, 80 pieces per pack
.021 x .025, 80 pieces per pack

Item-No. 620-1825N
Item-No. 620-2128N



Crimpable Stops or Spacers

Tubes .019 ID x .030 OD, 2 mm
40 pieces per pack
Tubes .030 ID x .041 OD, 2 mm
40 pieces per pack

Item-No. 620-0262
Item-No. 620-0263



Crimpable Hooks

Universal hooks for wires up to .018 x .025.

20 pieces per pack


Item-No. 620-3010



Hook Crimping Plier Angled

Designed for placement and securing of crimpable archwire hooks, stops and power hooks.

Item-No. 0L-204-333




Crimp Hooks TP – Slide Prove

Special interior coating effects a firm grip and prevents sliding.
 For wires up to .018 and .022 x .028.

10 pieces per pack

Item-No. 620-6010



Crimpable Stops TP – Slide Prove

Special interior coating effects a firm grip and prevents sliding.
 For wires up to .018 and .022 x .028.

10 pieces per pack

Item-No. 620-6011




Freddy Hook Buttons

NEW!

Smooth, low profile. Ideal for Lingual technology.

10 pieces per pack

Item-No. GE-FCB-001-01HK



Elast-O-Loop Ligatures


Package of 1008 rings.

Size .095, color: grey


7 stars à 144 rings per pack

Item-No. 630-6090

Latex-Free!



NEW



Standard

.090 EasyOn Safety Power Rings

- Flatter rectangular cross section .090 with soft edges.
- Tighter fit to very small lingual brackets.
- Increased elasticity and higher tear strength for a better control during a period of 6 – 8 weeks.
- Individual sticks reduce the risk of cross contamination.
- Resealable bags protect from harmful ultraviolet light and extend the shelf-life.

Silver: 40 sticks à 24 rings per pack

Item-No. 630-6510


Clear: 40 sticks à 24 rings per pack

Item-No. 630-6515

Black: 40 sticks à 24 rings per pack

Item-No. 630-6520

Separate leaflet for detailed information on request.




Elast-O-Loop Ligature .120

Latex free.

Black: 30 sticks à 24 rings per pack

Item-No. 630-0406



Elastomeric Ligature with Guard

Silver: 10 sticks à 10 rings per pack

Item-No. PP-406-429

Clear: 10 sticks à 10 rings per pack

Item-No. PP-406-430



Safety Ligature – Superstretch


Latex-Free!

Highly stretchable ligatures with a reset of more than 80%. Round rings, thin, particularly tear resistant and latex free.

Provides reduced friction with wires.

Black: 100 sticks à 10 rings per pack

Item-No. 630-6530



Memory Elastomeric Chain – Closed

Latex-Free!


Closed. All Japan Chains have patented formula which provides light continuous forces for weeks. Stain resistant and latex free.

Clear: 4.0 m per pack

Item-No. 640-0693

Grey: 4.0 m per pack

Item-No. 640-0698



Memory Elastomeric Chain – Narrow

Latex-Free!


Small rings, narrow space.

Clear: 4.7 m per pack

Item-No. 640-0691

Grey: 4.7 m per pack

Item-No. 640-0696



Memory Elastomeric Chain – Long

Latex-Free!


Large rings, long space.

Clear: 4.7 m per pack

Item-No. 640-0692

Grey: 4.7 m per pack

Item-No. 640-0697



Elastomeric Power Chain

Superior Elasticity and Rebound. Extremely High Resistance to Breakage and Tearing.

Clear Continuous 4.5 m per spool

Item-No. 400-05120

Clear Short 4.5 m per spool

Item-No. 200-05100

Clear Long 4.5 m per spool

Item-No. 300-05110

Grey Continuous 4.5 m per spool

Item-No. 400-11120

Grey Short 4.5 m per spool

Item-No. 200-11100

Grey Long 4.5 m per spool

Item-No. 300-11100

Black Continuous 4.5 m per spool

Item-No. 400-01120

Black Short 4.5 m per spool

Item-No. 200-01100

Black Long 4.5 m per spool

Item-No. 300-01110



Tooth Color or White Coated Ligatures – Short

White coated or tooth colored.
Ideal for ceramic brackets.


- .010 white 100 pieces per pack** [Item-No. 630-1210](#)
- .012 tooth color 100 pieces per pack** [Item-No. OL-100-412TC](#)
- .012 white 100 pieces per pack** [Item-No. OL-100-412TCW](#)



Preformed Ligature Wire – Long

Preformed Ligature Wire, long.

- .008 1000 pieces per pack** [Item-No. 630-1127](#)
- .009 1000 pieces per pack** [Item-No. 630-1130](#)
- .010 1000 pieces per pack** [Item-No. 630-1128](#)



Preformed Ligature Wire – Short

Preformed Ligature Wire, short.


- .008 1000 pieces per pack** [Item-No. 630-1140](#)
- .010 500 pieces per pack** [Item-No. 630-1145](#)



Kobayashi Twists – Short

Preformed Kobayashi Twists, short.

- .012 100 pieces per pack** [Item-No. OL-100-413](#)



Lingual Retainer Wire on Spool – Flat

Lingual Retainer Wire.
Flattened, 3-braided.

.010 x .022 spool with 1.5 m

Item-No. 650-1022

.012 x .025 spool with 1.5 m

Item-No. 650-1225



Gold Plated Retainer Wire – Multi Stranded

Stainless steel lingual retainer wire.
Gold coated 24 carat, 6-stranded.
Length: 355 mm



.0215 5 per pack

Item-No. 650-9020

.0195 5 per pack

Item-No. 650-9021

.0175 5 per pack

Item-No. 650-9023

.0155 5 per pack

Item-No. 650-9022



Bond-a-Braid Retainer Wire – Braided

Stainless steel lingual retainer wire.
Flattened, dead soft, 8-braided.
Length: 155 mm



.010 x .028 (0.25 x 0.711 mm)

Item-No. 650-8500

10 pieces per pack



Coaxial Wire – Straight

Stainless steel coaxial wire.
Length: 355 mm



.0150

Item-No. GH-MS6155

.0175

Item-No. GH-MS6175

.0195

Item-No. GH-MS6195

20 pieces per pack

Retainium™ Lingual Retainer – Single Strand


Nickel free single strand Titanium lingual retainer wire.
Length: 150 mm



.028 x .008 (0.711 x 0.211 mm)

Item-No. 650-8510

10 pieces per pack




Coaxial Wire on Spool

Stainless steel coaxial wire.

- .0150 spool with 30"**
- .0175 spool with 30"**
- .0195 spool with 30"**
- .0215 spool with 30"**

Item-No. GH-SST6155
Item-No. GH-SST6175
Item-No. GH-SST6195
Item-No. GH-SST6215



NiTi Superelastic on Spool

- .010 spool with 4.5 m**
- .012 spool with 3.0 m**
- .014 spool with 3.0 m**
- .016 spool with 3.0 m**

Item-No. 640-8977
Item-No. 640-8978
Item-No. 640-8979
Item-No. 640-8989


Stainless Steel Tubes – Round



Steel tube protects and prevents archwire from bending in case of large spaces. Protects the tongue from irritation. Length: 11.81" (300 mm) Inside diameter: 0.04" (1 mm)

1 piece per pack

Item-No. 620-4801



Protection Silicone Tube for Archs

Clear silicone tube for tongue protection. .018 inside diameter.

3 m per spool

Item-No. 642-0010

Cu-Alloy NiTi Archwires

Cu-Alloy wires consist of nickel, titanium and copper. The wires become active at mouth temperature and are flexible at room temperature. Smooth surface to reduce the friction between wires and brackets

10 pieces per pack

Wire Size	Ovoid Upper	Ovoid Lower	Nature Upper	Nature Lower	D-Form Univ.
.013	ORCU-213U	ORCU-213L	ORCU-313U	ORCU-313L	ORCU-913
.014	ORCU-214U	ORCU-214L	ORCU-314U	ORCU-314L	ORCU-914
.016	ORCU-216U	ORCU-216L	ORCU-316U	ORCU-316L	ORCU-916
.018	ORCU-218U	ORCU-218L	ORCU-318U	ORCU-318L	ORCU-918
.020	-	-	-	-	ORCU-920
.016 x .022	ORCU-21622U	ORCU-21622L	ORCU-31622U	ORCU-31622L	ORCU-91622
.014 x .025	ORCU-21425U	ORCU-21425L	ORCU-31425U	ORCU-31425L	ORCU-91425
.016 x .025	ORCU-21625U	ORCU-21625L	ORCU-31625U	ORCU-31625L	ORCU-91625
.017 x .025	ORCU-21725U	ORCU-21725L	ORCU-31725U	ORCU-31725L	ORCU-91725
.018 x .025	ORCU-21825U	ORCU-21825L	ORCU-31825U	ORCU-31825L	ORCU-91825
.019 x .025	ORCU-21925U	ORCU-21925L	ORCU-31925U	ORCU-31925L	ORCU-91925

NiTi Super Elastic Archwires

The nickel-titanium alloy allows specific properties, which make this wire the ideal treatment material. The super elastic NiTi wires operate at higher elasticity and lower stiffness with a very smooth and constant force effect. The excellent shape retention with the simultaneous resistance and predictable forces ensure fast results without the necessity of many adaptations or modifications. As a consequence, length of treatment is considerably reduced.

10 pieces per pack




Wire Size	Ovoid Upper	Ovoid Lower	Nature Upper	Nature Lower	D-Form Univ.
.012	ORNI-212U	ORNI-212L	ORNI-312U	ORNI-312L	ORNI-912
.014	ORNI-214U	ORNI-214L	ORNI-314U	ORNI-314L	ORNI-914
.016	ORNI-216U	ORNI-216L	ORNI-316U	ORNI-316L	ORNI-916
.018	ORNI-218U	ORNI-218L	ORNI-318U	ORNI-318L	ORNI-918
.020	ORNI-220U	ORNI-220L	ORNI-320U	ORNI-320L	ORNI-920
.016 x .016	ORNI-21616U	ORNI-21616L	ORNI-31616U	ORNI-31616L	ORNI-91616
.016 x .022	ORNI-21622U	ORNI-21622L	ORNI-31622U	ORNI-31622L	ORNI-91622
.017 x .022	ORNI-21722U	ORNI-21722L	ORNI-31722U	ORNI-31722L	ORNI-91722
.017 x .025	ORNI-21725U	ORNI-21725L	ORNI-31725U	ORNI-31725L	ORNI-91725
.018 x .022	ORNI-21822U	ORNI-21822L	ORNI-31822U	ORNI-31822L	ORNI-91822
.018 x .025	ORNI-21825U	ORNI-21825L	ORNI-31825U	ORNI-31825L	ORNI-91825
.019 x .025	ORNI-21925U	ORNI-21925L	ORNI-31925U	ORNI-31925L	ORNI-91925
.021 x .025	ORNI-22125U	ORNI-22125L	ORNI-32125U	ORNI-32125L	ORNI-92125

Stainless Steel Archwires

Manufactured in accordance with the highest standards of the industry. Strict production and quality controls ensure satisfaction with our products. For this archwire only medical 304V stainless steel is used.

10 pieces per pack

Wire Size	Ovoid Upper	Ovoid Lower	Nature Upper	Nature Lower	D-Form Univ.
.012	ORSS-212U	ORSS-212L	ORSS-312U	ORSS-312L	ORSS-912
.014	ORSS-214U	ORSS-214L	ORSS-314U	ORSS-314L	ORSS-914
.016	ORSS-216U	ORSS-216L	ORSS-316U	ORSS-316L	ORSS-916
.018	ORSS-218U	ORSS-218L	ORSS-318U	ORSS-318L	ORSS-918
.020	ORSS-220U	ORSS-220L	ORSS-320U	ORSS-320L	ORSS-920
.016 x .016	ORSS-21616U	ORSS-21616L	ORSS-31616U	ORSS-31616L	ORSS-91616
.016 x .022	ORSS-21622U	ORSS-21622L	ORSS-31622U	ORSS-31622L	ORSS-91622
.017 x .022	ORSS-21722U	ORSS-21722L	ORSS-31722U	ORSS-31722L	ORSS-91722
.017 x .025	ORSS-21725U	ORSS-21725L	ORSS-31725U	ORSS-31725L	ORSS-91725
.018 x .022	ORSS-21822U	ORSS-21822L	ORSS-31822U	ORSS-31822L	ORSS-91822
.018 x .025	ORSS-21825U	ORSS-21825L	ORSS-31825U	ORSS-31825L	ORSS-91825
.019 x .025	ORSS-21925U	ORSS-21925L	ORSS-31925U	ORSS-31925L	ORSS-91925
.021 x .025	ORSS-22125U	ORSS-22125L	ORSS-32125U	ORSS-32125L	ORSS-92125



NiTi Thermal Active Archwires

NiTi Thermal Active wires is Dead Soft when chilled below 5°C. The ideal force is achieved at 32°C. The wire distinguishes itself by the following properties: Smooth and consistent release of forces at a low and consistent force level as well as high deformation resistance.

10 pieces per pack

Wire Size	Ovoid Upper	Ovoid Lower	Nature Upper	Nature Lower	D-Form Univ.
.012	ORTN-212U	ORTN-212L	ORTN-312U	ORTN-312L	ORTN-912
.014	ORTN-214U	ORTN-214L	ORTN-314U	ORTN-314L	ORTN-914
.016	ORTN-216U	ORTN-216L	ORTN-316U	ORTN-316L	ORTN-916
.018	ORTN-218U	ORTN-218L	ORTN-318U	ORTN-318L	ORTN-918
.020	ORTN-220U	ORTN-220L	ORTN-320U	ORTN-320L	ORTN-920
.016 x .016	ORTN-21616U	ORTN-21616L	ORTN-31616U	ORTN-31616L	ORTN-91616
.016 x .022	ORTN-21622U	ORTN-21622L	ORTN-31622U	ORTN-31622L	ORTN-91622
.017 x .022	ORTN-21722U	ORTN-21722L	ORTN-31722U	ORTN-31722L	ORTN-91722
.017 x .025	ORTN-21725U	ORTN-21725L	ORTN-31725U	ORTN-31725L	ORTN-91725
.018 x .022	ORTN-21822U	ORTN-21822L	ORTN-31822U	ORTN-31822L	ORTN-91822
.018 x .025	ORTN-21825U	ORTN-21825L	ORTN-31825U	ORTN-31825L	ORTN-91825
.019 x .025	ORTN-21925U	ORTN-21925L	ORTN-31925U	ORTN-31925L	ORTN-91925
.021 x .025	ORTN-22125U	ORTN-22125L	ORTN-32125U	ORTN-32125L	ORTN-92125

Beta Titanium Archwires

The quality of Beta Titanium wires lies between the bending properties of stainless steel wires and the flexibility of NiTi wires. Its hardness is about 42% of stainless steel wires and is easy to bend and easy to handle. It is perfectly suitable in the final treatment phase to apply stops or loops.

10 pieces per pack


Wire Size	Ovoid Upper	Ovoid Lower	Nature Upper	Nature Lower
.016	ORBT-216U	ORBT-216L	ORBT-316U	ORBT-316L
.018	ORBT-218U	ORBT-218L	ORBT-318U	ORBT-318L
.020	ORBT-220U	ORBT-220L	ORBT-320U	ORBT-320L
.016 x .016	ORBT-21616U	ORBT-21616L	ORBT-31616U	ORBT-31616L
.016 x .022	ORBT-21622U	ORBT-21622L	ORBT-31622U	ORBT-31622L
.017 x .022	ORBT-21722U	ORBT-21722L	ORBT-31722U	ORBT-31722L
.017 x .025	ORBT-21725U	ORBT-21725L	ORBT-31725U	ORBT-31725L
.018 x .018	ORBT-21818U	ORBT-21818L	ORBT-31818U	ORBT-31818L
.018 x .025	ORBT-21825U	ORBT-21825L	ORBT-31825U	ORBT-31825L
.019 x .025	ORBT-21925U	ORBT-21925L	ORBT-31925U	ORBT-31925L
.021 x .025	ORBT-22125U	ORBT-22125L	ORBT-32125U	ORBT-32125L


RCS NiTi Superelastic Archwires

RCS NiTi Superelastic offers superb rebound and memory characteristics to ensure desired results without the need for adjustment or modification. It allows a high degree of deflection without causing deformation of the original arch form characteristics. Efficiency of performance ensures achievement of objectives without the need for multiple staging of archwires as with stainless steel alternatives.

10 pieces per pack

Wire Size	RCS Upper	RCS Lower
.012	ORRN-412U	ORRN-412L
.014	ORRN-414U	ORRN-414L
.016	ORRN-416U	ORRN-416L
.018	ORRN-418U	ORRN-418L
.020	ORRN-420U	ORRN-420L
.016 x .016	ORRN-41616U	ORRN-41616L
.016 x .022	ORRN-41622U	ORRN-41622L
.017 x .022	ORRN-41722U	ORRN-41722L
.017 x .025	ORRN-41725U	ORRN-41725L
.018 x .022	ORRN-41822U	ORRN-41822L
.018 x .025	ORRN-41825U	ORRN-41825L
.019 x .025	ORRN-41925U	ORRN-41925L
.021 x .025	ORRN-42125U	ORRN-42125L





Dry Field System Standard

Absolutely necessary in lingual technique.
Autoclavable, not exceed 280°F (135°C).
Color: white

1 piece per pack

Item-No. 735-0010



Dry Field System Small

Absolutely necessary in lingual technique.
Autoclavable, not exceed 280°F (135°C).
Color: red

1 piece per pack

Item-No. 735-0009



Dry Field Tongue Guards

Replacement tongue guards for Dry Field System.
Clear: large; Red: small

Large: 2 pieces per pack

Item-No. 735-0011

Small: 2 pieces per pack

Item-No. 735-0012



Dry Field Low Volume Adapters

Y-Connector for the Dry Field System.

4 pieces per pack

Item-No. 735-0006



Dry Field Replacement Silicon Tubes

Silicon replacement tubes for the Dry Field System.

2 pieces per pack

Item-No. 735-0015

10 pieces per pack

Item-No. 735-0015/10



Dry Tips

Saliva absorbent with impermeable plastic film.
For use with Dry Field System.

Large: 50 pieces per pack

Item-No. 735-0051

Small: 50 pieces per pack

Item-No. 735-0050



Photo Cheek Retractor Double Ended

Clear. Set small and medium. For the posterior or lateral area.
Can only be cold sterilized.

2 pieces per pack

Item-No. 720-0006



Photo Cheek Retractor Double Ended – Soft

Opaque. Set small and medium. Autoclavable up to 121°C

2 pieces per pack

Item-No. 720-0006S



Photo Cheek Retractor Double Ended

Clear. Set medium and large. For the posterior or lateral area.
Can only be cold sterilized.

2 pieces per pack

Item-No. 720-0007



Photo Cheek Retractor Double Ended – Soft

Opaque. Set medium and large. Autoclavable up to 121°C

2 pieces per pack

Item-No. 720-0007S



Photo Cheek Retractor Double Ended

Autoclavable up to 121°C.

2 pieces per pack

Item-No. 0L-300-027



Photo Cheek Retractor Double Ended – Small

Clear. Cold sterilizable
and autoclavable up to 121°C.

2 pieces per pack

Item-No. 720-0008




Photo Cheek Retractor Single Side


Transparent Blue. Cold sterilizable.

**Lip Retractor,
2 pieces per pack**

Item-No. HW-605450

**Mini Lip Retractor,
2 pieces per pack**

Item-No. HW-605451




Lip Retractor

Clear. For occlusal photographs.
Cold sterilizable or autoclavable up to 121°C.

Small: 1 piece per pack
Large: 1 piece per pack

Item-No. 720-0012
Item-No. 720-0011




Lip Retractor Long Handle

Clear. Cold sterilizable
and autoclavable up to 121°C.

Small: 1 piece per pack
Large: 1 piece per pack

Item-No. 720-0010
Item-No. 720-0009




Lip Retractor Long Handle – Soft

Opaque. Cold sterilizable
and autoclavable up to 121°C.

Small: 1 piece per pack
Large: 1 piece per pack

Item-No. 720-0010S
Item-No. 720-0009S



Cheek Retractor Small


Polycarbonate. Only cold sterilization.

 **Red**

Item-No. 730-0002

1 piece per pack

Other sizes and models available.




Cheek Retractor Medium – Soft

Polypropylene. Autoclavable up to 121°C.

 **Red**

Item-No. 730-0010S


 **Opaque**

Item-No. 730-0013S

 **Blue**

Item-No. 730-0015S

1 piece per pack



Cheek Retractor Medium

Polycarbonate. Only cold sterilization.

 **Red**

Item-No. 730-0010

 **White**

Item-No. 730-0011

 **Clear**

Item-No. 730-0013

 **Blue**


Item-No. 730-0015

 **Yellow**

Item-No. 730-0017

1 piece per pack

Other sizes and models available.



Cheek Retractor Large

Polycarbonate. Only cold sterilization.

 **Blue**

Item-No. 730-0025


 **Red**

Item-No. 730-0026

 **Yellow**

Item-No. 730-0027

1 piece per pack




Cheek Retractor Small

Clear. Only cold sterilization.

2 pieces per pack

Item-No. OL-300-021




Cheek Retractor Medium

Clear. Only cold sterilization.

2 pieces per pack

Item-No. OL-300-022



Cheek Retractor Large

Clear. Only cold sterilization.

2 pieces per pack

Item-No. OL-300-020



Cheek Retractor Assortment

Cheek Retractor Assortement with 3 sizes (large, medium and small).

3 pieces per pack (1 per size)

Item-No. OL-300-023



Orolin Perasept

Concentrated powder for cold sterilization of heat-sensitive instruments and cheek retractors.

1 kg per pack

Important!

Item-No. 710-0500




Cheek Retractor Stainless Steel Double Ended

Autoclavable, chemiclavable
and sterilisable with hot air.

2 pieces per pack

Item-No. 300-16-216B



Cheek Retractor Steel Sternberg

Autoclavable, chemiclavable
and sterilisable with hot air.

2 pieces per pack

Item-No. 730-0060

Rhodium Coated Mirrors




Photo Mirror with Handle – Small

Photo mirrors with functional handle.
For occlusal photographs. Size: 75 x 65 mm.
With rhodium coated glass.

1 piece per pack

Item-No. 720-0017




Photo Mirror with Handle – Medium

Photo mirrors with functional handle.
For occlusal photographs. Size: 90 x 70 mm.
With rhodium coated glass.

1 piece per pack

Item-No. 720-0018




Photo Mirror with Handle – Large

Photo mirrors with functional handle.
For occlusal photographs. Size: 105 x 75 mm.
With rhodium coated glass.

1 piece per pack

Item-No. 720-0019





Photo Mirror with Handle – Lateral Photos

Photo mirrors with functional handle.
For lateral photographs. Size: 125 x 50 mm.
With rhodium coated glass.

1 piece per pack

Item-No. 720-0020

**Please pay attention to
the sterilization instructions!**



Ultra FS High Tech Mirrors

High Tech for High Quality

Superior reflection in terms of sharpness, natural colors and brightness. Developed and optimized by ion-beam-curing, our Photo Mirror achieve unsurpassed reflection and coating hardness.

- Exactly sharp, distortion-free, high-contrast image to the smallest detail.
- Brilliant color naturally and pleasantly light at an optimal uniform reflectance spectrum.
- Highly resistant to all kinds of disinfection and sterilization (180°C) to acidic refresh cleaners.




Photo Mirror Occlusal Ultra FS Surface – Small

Size: 76 x 58 mm

1 piece per pack

Item-No. 721-0017




Photo Mirror Occlusal Ultra FS Surface – Medium

Size: 95 x 66 mm

1 piece per pack

Item-No. 721-0018




Photo Mirror Occlusal Ultra FS Surface – Large

Size: 106 x 74 mm

1 piece per pack

Item-No. 721-0019




Photo Mirror Lateral Ultra FS Surface – Small

Size: 120 x 36 mm

1 piece per pack

Item-No. 721-0020




Photo Mirror Lateral Ultra FS Surface – Large


Size: 128 x 42 mm

1 piece per pack

Item-No. 721-0021



Chrome Plated Mirrors




Chrome Plated Mirror – Palatal Adult

Photo mirrors for sharp and clear images.
For occlusal photographs. Size: 59/69 x 103 mm.
With chrome plated surface.

1 piece per pack

Item-No. OL-300-001




Chrome Plated Mirror – Palatal Child

Photo mirrors for sharp and clear images.
For occlusal photographs. Size: 51/61 x 88 mm.
With chrome plated surface.

1 piece per pack

Item-No. OL-300-002




Chrome Plated Mirror – Lingual Narrow

Photo mirrors for sharp and clear images.
For occlusal photographs. Size: 36 x 151 mm.
With chrome plated surface.

1 piece per pack

Item-No. OL-300-003




Chrome Plated Mirror – Lingual Wide

Photo mirrors for sharp and clear images.
For occlusal photographs. Size: 39/45 x 154 mm.
With chrome plated surface.

1 piece per pack

Item-No. OL-300-004



Handle for Mirrors

Handle for photo mirrors up to 3 mm thickness.

1 piece per pack

Item-No. OL-GH-SGPR



Plicafol

Adaptable space maintaining foil.

100 pieces per pack

Item-No. 830-3910



Bisico S1 suhy*

Kneadable preliminary impression material for the first step in the classic double impression technique.

2 x 300 cc per pack

Item-No. 830-1090



Bisico Adhesive*

This tray adhesive improves the reliable adhesion of the impression material to the tray.

10 g bottle

Item-No. 830-8210



Bisico S4 suhy*

Addition-curing, extremely fine-flowing correction impression material with superhydrophilic properties for use in all impression techniques.

**3 x 50 cc, 18 mixing tips,
10 intraoral tips per pack**

Item-No. 830-7143



Bisico Dispenser 50*

For the new cartridge system.

1 piece per pack

Item-No. 830-1910



Bisico Mixing Tips*

Mixing tips, yellow.

50 pieces per pack

Item-No. 830-1950



Bisico IOT Intraoral Mixing Tips*

For mixer, yellow.

50 pieces per pack

Item-No. 830-1960

*Not sold in France

Improved formula

Dreve Dynax® putty



Dental impression material on the basis of vinylpolysiloxane.

Advantages

- Non-sticky, easily kneadable consistency
- Suitable for sandwich-, correction-, double mixing technique and valve border impressions
- Especially suitable for patients with remaining teeth suffering from periodontitis
- Scannable surface
- High hydrophilic properties
- Innovative snap-effect
- Outmatched mechanical properties

Technical Data

Colour	sunny yellow
Final hardness	60 Shore A
Consistency	kneadable type 0
Processing time	90 sec.
Removal	after 90 sec.
Mixing ratio of components A+B	1:1

2 x 450 ml (1 x A/1 x B) incl. dosing spoons [Item-No. D6101](#)

Dreve Dynax® light



Addition-curing, thin flowing correction silicone with outstanding detail precision and perfect thixotropy for the correction technique.

2 x 50 ml Double Cartridge

[Item-No. D6103](#)

8 x 50 ml Double Cartridge

[Item-No. D6104](#)

Dreve Dynax® mono



Addition-curing precision impression silicone.
A perfect impression-taking with only one material!

2 x 50 ml Double Cartridge

[Item-No. D6106](#)

8 x 50 ml Double Cartridge

[Item-No. D6107](#)



Dreve Dynax® heavy body

Addition-curing „tray type“ silicone.
A particularity in our range of products and
a real alternative to polyether impressions.

2 x 50 ml Double Cartridge
8 x 50 ml Double Cartridge

Item-No. D6112
Item-No. D6113



Dreve Dynax® clear

A real innovation worldwide: Clear, transparent,
addition-curing, monophasic impression silicone.
Ensures a perfect control already during the
impression-taking. Universally applicable and
particularly suitable for the implant impression-taking.

2 x 50 ml Double Cartridge
8 x 50 ml Double Cartridge

Item-No. D6121
Item-No. D6122



Dreve Injector

Dosing unit for Dynax® light double
cartridges with mixing ratio 1:1.

Injector white: 1 piece per pack

Item-No. D5007



Dreve Mixing Tips

Mixing tip for Dynax® light.

- Mixing ratio 1:1
- Inner diameter 4.2 mm
- Color: yellow

40 pieces per pack

Item-No. D4970X5



Dreve Application Tips

For mixing tips and intra-oral application.

40 pieces per pack

Item-No. D49604



Terry Protection Towel


Patient face-mask. Covers the face during the work with the Microetcher. Washable at 40°C.

White: 12 pieces per pack

Item-No. 735-0060

Light Blue: 12 pieces per pack

Item-No. 735-0061



Monoart Capri Towel Blue

To cover the patient's face. One-sided film coating. Size: 33 x 22 x 8 cm

500 pieces per pack

Item-No. 735-0500



Mirafloss

Very tear-resistant.

Big Box Dispenser with 2 x 200 m spools.

White, waxed

Item-No. 811-5735

Refill:

1 x 200 m spool;

White, waxed

Item-No. 811-5736




Occlusions Foil

Occlusions foil BK 05, Bausch KG.

200 foils per pack

Item-No. 735-0000

DuraLock®



DuraLock®

- | | |
|------------------------|-------------------------|
| 1: Upper red | Item-No. 640-041 |
| 2: Upper orange | Item-No. 640-042 |
| 3: Upper blue | Item-No. 640-043 |
| 4: Upper green | Item-No. 640-044 |
| 5: Upper yellow | Item-No. 640-045 |
| 6: Upper purple | Item-No. 640-046 |

- | | |
|------------------------|-------------------------|
| 1: Lower red | Item-No. 640-051 |
| 2: Lower orange | Item-No. 640-052 |
| 3: Lower blue | Item-No. 640-053 |
| 4: Lower green | Item-No. 640-054 |
| 5: Lower yellow | Item-No. 640-055 |
| 6: Lower purple | Item-No. 640-056 |

10 pieces per pack

DuraLock® Starter Kit

Assorted colors **Item-No. 640-040/1**

12 pieces per pack (1 each size)

Only Cold sterilization.
For sterilization we recommend
Orolin Perasept on page E5.



Retainer Box

Colorful mixed with glitter
(orange, black, green, pink, blue).
Inner height: 2.4 cm.

10 pieces per pack

Item-No. OL-700-0006



Brace Gard

Silicon protection for lips and tongue.
Mini mirror jars. Box with 3 strips.

30 boxes per pack

Item-No. 735-0037




Utility Wax Bars

Wax bars in boxes.

435 g per pack

Item-No. 820-9010



Happy Morning

Single toothbrush with paste.

100 pieces per pack

Item-No. OL-999-0101




Chromat Alginate

Alginate impression material.
 Setting time 2 min.
 Low dust, lead and cadmium free.
 Excellent biocompatibility.
 Chromatic phase indicator.

Scale Prices

1 x 450 g pouch per pack
24 x 450 g pouch per pack
48 x 450 g pouch per pack










Item-No. 586-0011



Model Storage Boxes without frames for Labels

These boxes have been specifically designed to store orthodontic models. Stackable. Available in 9 different colors. Model without frames for labels. Size: 305 x 85 x 67 mm.

25 pieces per pack










	White	589-8010		Brown	589-8015
	Blue	589-8011		Grey	589-8016
	Red	589-8012		Beige	589-8017
	Green	589-8013		Ivory	589-8018
	Yellow	589-8014			




Model Storage Boxes with frames for Labels

These boxes have been specifically designed to store orthodontic models. Stackable. Available in 9 different colors. Model with frames for labels. Size: 305 x 85 x 67 mm.

25 pieces per pack / Labels: 25 pieces per pack

	White	589-8010E		Brown	589-8015E
	Blue	589-8011E		Grey	589-8016E
	Red	589-8012E		Beige	589-8017E
	Green	589-8013E		Ivory	589-8018E
	Yellow	589-8014E		Labels	589-8020



Model Storage Boxes – Tackable or Foldable

Stable, stackable, easy labeling. Available in two versions: foldable and tackable. Shipped in single flat pieces for convenient storage. Size: 315 x 93 x 70 mm.

Scale Prices

A) Tackable: 100, 500 or 1000 pieces **Item-No. 100-0000**
B) Foldable: 100, 500 or 1000 pieces **Item-No. 101-0000**

The Aqualizer™ - The Original

Water splint relaxes muscles and relieves pain


The Aqualizer™ is the only hydrostatic bite-splint, which can be immediately incorporated to each patient without preparation or laboratory work. The Aqualizer™'s built-in fluid-system adjusts automatically to the patient's anatomy. It allows you to treat your patient's pain or stiffness quickly and accurately. For primary care and introduction to a therapie, the Aqualizer™ has proven effective.

The Aqualizer™ offers three fluid volumes: Low (1 mm), Medium (2 mm) and High (3 mm). The volume of fluid in the Aqualizer™ controls the vertical dimension. All fluid volumes cover the same arch length.



The Aqualizer™ supports the treatment of typical TMD syndrome symptoms

- Persistent TMJ pain and headaches. The pain can be mild or severe and can start in the front of the head, along the sides, in the upper part of the neck, or it can radiate to all areas.
- Ear pain or stuffiness and tinnitus.
- Facial pain or tightness.
- Limited ability to open the mouth.
- Pain above or behind the eyes, neck and shoulder pain



For more information ask for our detailed brochure.

Aqualizer™ Slim (Classic)

The original design. Can cause irritation in susceptible patients at the gingival but is resistant to strong biting forces.

High, 3 mm fluid volume **Item-No. Z-AQ-300**

Medium, 2 mm fluid volume **Item-No. Z-AQ-301**

Low, 1 mm fluid volume **Item-No. Z-AQ-304**

Aqualizer™ Ultra

The Aqualizer™ Ultra is a new improved version of the Aqualizer™ designed to be more comfortable to the gums than previous models. This is our bestseller and will fit most adult patients' mouths.

High, 3 mm fluid volume **Item-No. Z-AQ-305**

Medium, 2 mm fluid volume **Item-No. Z-AQ-306**

Low, 1 mm fluid volume **Item-No. Z-AQ-307**

Aqualizer™ Mini

High, 3 mm fluid volume **Item-No. Z-AQ-308**

Medium, 2 mm fluid volume **Item-No. Z-AQ-309**

Low, 1 mm fluid volume **Item-No. Z-AQ-310**

Ordering Information



Carmen Moscato

Order Processing, Customer Service

Fon: +49 (0)8381/9485937

Fax: +49 (0)8381/9485939

E-Mail: c.moscato@pelzundpartner.net



Ingrid Fischer

Order Processing, Customer Service

Fon: +49 (0)8381/8909524

Fax: +49 (0)8381/9485939

E-Mail: ingrid.fischer@pelzundpartner.net



Elke Buhmann

Order Processing, Customer Service

Fon: +49 (0)8381/8909527

Fax: +49 (0)8381/9485939

E-Mail: e.buhmann@pelzundpartner.net



Elsa Cid Berger

Order Processing, Customer Service

Fon: +49 (0)8381/8909528

Fax: +49 (0)8381/9485939

E-Mail: e.berger@pelzundpartner.net



Hubert Stiefenhofer

Managing Director

Fon: +49 (0)8381/8909511

E-Mail: h.stiefenhofer@pelzundpartner.net



Brigitte Koch

Accounts Department (**Mo – Th, 1 to 5 p.m.**)

Fon: +49 (0)8381/8909534

Fax: +49 (0)8381/9485939

E-Mail: b.koch@pelzundpartner.net

We agree to 5% discount for orders superior to 250 Euro net values,
and to 10% discount for orders superior to 500 Euro (Excluded Online-Shop!).

**Do you need help or have questions
about our products?**

**Call us – we will be glad to help
and advise you.**

Your Pelz & Partner team

Since 2008

Pelz & Partner[®]
Specialist for Lingual Technique and Orthodontics



Pelz & Partner[®]

Specialist for Lingual Technique and Orthodontics



Pelz & Companion GmbH
Lauenbuehlstrasse 59
D-88161 Lindenberg/Allg.
Germany

Fon +49 (0) 83 81 948 59 37

Fax +49 (0) 83 81 948 59 39

E-Mail info@pelzundpartner.net

Web www.pelzundpartner.net

